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
### ERA (ENDOMETRIAL RECEPTIVITY ANALYSIS)

Patient information	Sample information	Clinic information
Unique pat id.:	Date received:	Clinic:
Sample type:	Report Date:	Clinician: Dr.
Patient name:	First intake of P4:	No. biopsy:
Patient DOB:	Date of biopsy:	
	Cycle type:	

**TEST RESULTS:**

**PRE-RECEPTIVE**

**Recommendation: The personalized embryo transfer (pET) of a blastocyst/s should be performed with 146 ± 3 hours of progesterone administration (1 day later than the time at which this endometrial biopsy was performed). A new endometrial biopsy is not required. \*\***

**INTERPRETATION OF YOUR RESULT:**

According to our internal data, 89% of women with similar endometrial profile reached receptivity with 1 more day of progesterone administration (confidence interval of 95% [86%-91%]), so in these cases new endometrial biopsy is not needed. Therefore, blastocyst/s transfer is recommended with 146 ± 3 hours of progesterone administration.

For a day-3 embryo/s, the transfer should be performed two days earlier than indicated in the recommendation for blastocyst transfer above.

\*\* This recommendation is only applicable to the same type of cycle treatment as the one used for this endometrial biopsy and if the endogenous progesterone measured prior to the first progesterone intake is <1ng/ml.


**TEST DESCRIPTION:**

ERA (Endometrial Receptivity Analysis) is a molecular tool used to determine if the endometrium (the mucous membrane lining the womb) exhibits a receptive profile after 5 days of progesterone exposure, the time at which the endometrium is typically ready for embryo implantation. This molecular diagnosis method is based on measuring the gene expression profile of endometrial tissue. Therefore, ERA helps to determine when the endometrium presents the ideal condition for embryo implantation, increasing the possibility of a successful in vitro fertilization treatment.

**COMMENTS**

None

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## Referências

Díaz-Gimeno P, Horcajadas JA, Martínez-Conejero JA, Esteban FJ, Alamá P, Pellicer A, Simón C. A genomic diagnostic tool for human endometrial Receptivity based on the transcriptomic signature. *Fertil Steril.* 2011; 95(1):50-60, 60.e1-15.

Díaz-Gimeno P, Ruiz-Alonso M, Blesa D, Bosch N, Martínez-Conejero JA, Alamá P, Garrido N, Pellicer A, Simón C. The accuracy and reproducibility of the endometrial Receptivity array is superior to histology as a diagnostic method for endometrial Receptivity. *Fertil Steril.* 2013; 99(2):508-17.

Clemente-Ciscar M, Ruiz-Alonso M, Blesa D, Jimenez-Almazan J, Bahceci M, Banker M et al. Endometrial receptivity analysis (ERA) using a next generation sequencing (NGS) predictor improves reproductive outcome in recurrent implantation failure (RIF) patients when compared to ERA arrays. *Hum Reprod.* 2018; 33(Supp1):8-8.

Findikli N, Gultomruk M, Boynukalin K, Kavrut M, Oral E, Bahceci M. Combinatorial use of Endometrial Receptivity Array (ERA) and PGT-A can improve the clinical outcome in cases with previous ART failures. *Hum Reprod.* 2018; 33(Supp1):84-85.

# ERA<sup>®</sup>

Endometrial  
Receptivity Analysis

Hashimoto T, Koizumi M, Doshida M, Toya M, Sagara E, Oka N. et al. Efficacy of the endometrial Receptivity Array for repeated implantation failure in Japan: A retrospective, two-centers study. *Reprod Med Biol.* 2017; 16(3): 290-296.

Mahajan N. Endometrial receptivity array: Clinical application. *J Hum Reprod Sci.* 2015; 8(3):121-9.

Pasternak M, Schattman G, Rosenwaks Z. Pregnancy outcomes in patients undergoing embryo transfer in cycle following endometrial Receptivity assay. *Fertil Steril.* 2018; 110(4):e243-244.

Ruiz-Alonso M, Blesa D, Díaz-Gimeno P, Gómez E, Fernández-Sánchez M, Carranza F, Carrera J, Vilella F, Pellicer A, Simón C. The endometrial Receptivity array for diagnosis and personalized embryo transfer as a treatment for patients with repeated implantation failure. *Fertil Steril.* 2013; 100(3):818-24.

Ruiz-Alonso M, Galindo N, Pellicer A, Simón C. What a difference two days make: “personalized” embryo transfer (pET) paradigm: a case report and pilot study. *Hum Reprod.* 2014; 29(6):1244-7.

Taguchi S, Funabiki M, Hayashi T, Tada Y, Iwaki Y, Karita M, et al. The implantation rate of Japanese infertile patients with repeated implantation failure can be improved by endometrial Receptivity array (ERA) test: A randomized controlled trial. *Fertil Steril.* 2018; 110(4):e90.

# EMMA

Endometrial Microbiome  
Metagenomic Analysis

# EMMA

Endometrial Microbiome  
Metagenomic Analysis

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EMMA TEST RESULT ABNORMAL ENDOMETRIAL MICROBIOME.  
ALICE TEST RESULT NEGATIVE FOR BACTERIAL PATHOGENS CAUSING CHRONIC ENDOMETRITIS

EMMA		ALICE	
Most abundant bacteria	%	Chronic Endometritis pathogens	%
<i>Lactobacillus</i>	22.86% *	<i>Enterobacteriaceae</i>	Not Detected
<i>Gardnerella</i>	36.43%	<i>Escherichia</i>	Not Detected
<i>Prevotella</i>	22.55%	<i>Klebsiella</i>	Not Detected
<i>Bifidobacterium</i>	13.64%	<i>Enterococcus</i>	Not Detected
Others	4.52%	<i>Chlamydia</i>	Not Detected
		<i>Mycoplasma</i>	Not Detected
		<i>Neisseria</i>	Not Detected
		<i>Ureaplasma</i>	Not Detected
		<i>Streptococcus</i>	Not Detected
		<i>Staphylococcus</i>	Not Detected

\* For reference intervals, please refer to Moreno et al. Am J Obstet Gynecol. 2018.

Color	Genus
Blue	Gardnerella
Orange	Lactobacillus
Dark Blue	Prevotella
Light Blue	Bifidobacterium
Grey	Others

#### INTERPRETATION OF YOUR RESULT AND RECOMMENDATION

DNA from bacterial pathogens of the reproductive tract has been detected in a significant amount in the endometrial sample.

Antibiotic therapy followed by probiotic treatment is recommended before continuing with ART. Please find below the suggested therapy based on the bacteria detected. We also recommend the analysis of a second sample after treatment, to confirm the restoration of a favourable environment for implantation.

#### SUGGESTED THERAPY


Metronidazole 500mg/12h for 7 days followed by probiotic treatment is recommended. A list with recommended probiotics of vaginal administration is provided at the end of this report.

Support +34 916993100 IGENOMIX SPAIN LAB, S.L.U. C/ Nave de Masamiron 11, Distrito Barajas B, Parque Tecnológico, ES-28960 PATERNA, Valencia, Tef: +34 916 995 310 MEE-0026024 - 115  
Igenomix is present in: Los Angeles | Miami | Sao Paulo | Valencia | Dubai | New Delhi | Mexico | New York | Montreal | London | Tokyo | Buenos Aires | Marostica | Istanbul SPA\_3\_2\_25\_254\_001-25

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## Referências

Human Microbiome Project Consortium. Structure, function and diversity of the healthy human microbiome. *Nature* 2012; 486:207–14.

Moreno I, Codoñer FM, Vilella F, Valbuena D, Martinez-Blanch JF, Jimenez-Almazan J, Alonso R, Alama P, Remohi J, Pellicer A, Ramon D, Simon C. Evidence that the endometrial microbiota has an effect on implantation success or failure. *Am J Obstet Gynecol.* 2016; 215:684-703.

Inmaculada Moreno; Jason M Franasiak. Endometrial microbiota - new player in town. *Fertility and Sterility.* 2017;108, pp. 32 - 39.

Jason M Franasiak; Inmaculada Moreno; Carlos Simon. Microbiome in Embryonic Implantation and Implantation Failure. In: *Recurrent Implantation Failure, Etiologies and Clinical Management.* 2018; Chapter 11, pp. 175 - 195. Springer, Cham. ISBN 978-3-319-71966-5

Carlos Simon; Inmaculada Moreno. Deciphering the effect of reproductive tract microbiota on human reproduction. *Reproductive Medicine and Biology.* 2018;18 - 1, pp. 40 - 50.

Inmaculada Moreno; Carlos Simon. Relevance of assessing the uterine microbiota in infertility. *Fertility and Sterility.* 2018 Aug;110(3):337-343. doi: 10.1016/j.fertnstert.2018.04.041.

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### ANALYSIS OF INFECTIOUS CHRONIC ENDOMETRITIS (ALICE)

Patient information	Sample information	Clinic information
Unique pat id:	Date received:	Clinic:
Patient name:	Report date/time:	Clinician:
Patient DOB:	Sample type: Endometrial Biopsy	
Allergic to antibiotics	Cycle type:	
	No. Biopsy:	
	Date of biopsy:	

ALICE TEST RESULT

POSITIVE FOR BACTERIAL PATHOGENS CAUSING CHRONIC ENDOMETRITIS

Chronic Endometritis pathogens	%
<i>Enterobacteriaceae</i>	Not Detected
<i>Escherichia</i>	Not Detected
<i>Klebsiella</i>	Not Detected
<i>Enterococcus</i>	Not Detected
<i>Chlamydia</i>	Not Detected
<i>Mycoplasma</i>	Not Detected
<i>Neisseria</i>	Not Detected
<i>Ureaplasma</i>	Not Detected
<i>Staphylococcus</i>	Not Detected
<i>Streptococcus</i>	69,57%

### INTERPRETATION OF YOUR RESULT AND RECOMMENDATION

DNA from bacterial pathogens of the reproductive tract causing chronic endometritis has been detected in a significant amount in the endometrial sample.

In this case, the same pathogens detected in the previous ALICE test are still present in the sample after treatment. This could be due to various reasons including an incomplete therapy or resistance to the administered antibiotic. For this reason, a second line therapy with a different antibiotic followed by probiotic is recommended below. We also recommend the analysis of another biopsy after treatment to confirm the restoration of a physiological environment.

### SUGGESTED THERAPY


A second line treatment recommended would be Diacetyl-Midecamycin 600mg/12h for 7 days by oral way followed by vaginal probiotic treatment. A list with recommended probiotics of vaginal administration is provided at the end of this report.

# ALICE

Analysis of Infectious  
Chronic Endometritis

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Cicinelli E, Matteo M, Tinelli R, Pinto V, Marinaccio M, Indraccolo U, De Ziegler D, Resta L. Chronic endometritis due to common bacteria is prevalent in women with recurrent miscarriage as confirmed by improved pregnancy outcome after antibiotic treatment. *Reprod Sci* 2014; 21(5):640-7.

Cicinelli E, Matteo M, Tinelli R, Lepera A, Alfonso R, Indraccolo U, Marocchella S, Greco P, Resta L. Prevalence of chronic endometritis in repeated unexplained implantation failure and the IVF success rate after antibiotic therapy. *Hum Reprod*, 2015; 30(2):323-30.

# ALICE

Analysis of Infectious  
Chronic Endometritis

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